

Anti-HA (H₁N₁/California/06/2009)

Product Name	Anti-HA (H ₁ N ₁ /California/06/2009)
Catalog No	IT-003-SW
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HA (H ₁ N ₁ /California/06/2009) protein (Genebank accession ACP41935)
Applications	Western Blotting (1:200-1:500)
Specificity	Reacts with HA (H ₁ N ₁ /California/06/2009) protein. Cross-reactivity to other subtypes not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

