



# Anti- H1(A/swine/Pan)

Product Name	Anti- H1(A/swine/Pan)
Catalog No	IT-003-SW-Pan-M1
Source	Mouse monoclonal antibody, clone 43B10
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing H1(A/swine/Guangxi/13/2006(H1N2)) (Genebank No. <a href="#">ABQ42444</a> )
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with most H1 of swine origin. Not cross-react with other H1 and other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-003-SW-Pan-M1

