




Anti-HBeAg (HBV)

Product Name	Anti-HBeAg (HBV)
Catalog No	IT-004-003
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Hepatitis B virus HBeAg antigen (Genebank accession #AAO63518)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HBV HBeAg proteins. Cross-reactivity to other subtypes not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control;
E: 293 cell expressing HBeAg antigen (HBV)

Note: