

# Anti-E2(HCV)



Product Name	Anti-E2(HCV)
Catalog No	IT-004-0052
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Peptide “GPVYCFTPSPVVVGTTD” from HCV E2 protein (a.a. 119-135). This epitope is generally conserved in all genotypes.
Applications	ELISA(1:1000-1:100,000); WB(1:200-1:2000), IP, etc.
Specificity	Reacts with E2 protein from most genotypes of HCV.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

