



# Anti-NS3(HCV)

Product Name	Anti-NS3 (HCV)
Catalog No	IT-004-011M7
Source	Mouse monoclonal antibody
Isotype	IgM, clone 10B4
Formulation	Each vial contains 100 µg purified IgM in 100 µl PBS with 40% glycerol
Immunogen	Peptide from a conserved region of NS3 polypeptide of Hepatitis C virus (a.a. 1458-1467) (CNTCVTQTVD)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with a conserved epitope in NS3 polypeptide of Hepatitis C virus (a.a. 1458-1467) (CNTCVTQTVD)
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-004-011M7**

