

Anti-NS3 (HCV)



Product Name	Anti- NS3 (HCV)
Catalog No	IT-004-014
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Peptide from a conserved region of NS3 polypeptide of Hepatitis C virus (a.a. 1617-1628) (LIRLKPTLHGPT)
Applications	WB(1:2000), ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with NS3 polypeptide of Hepatitis C virus (a.a. 1617-1628) (LIRLKPTLHGPT).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-004-014

