

Anti-HDAg(2a)



Product Name	Anti- HDAg(2a)
Catalog No	IT-004-021EM3
Source	Mouse monoclonal antibody, clone 1G11
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant small delta antigen (Hepatitis delta virus) protein (IT-004-021Ep).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HDAg(2a). Cross-reactivity not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

