

# P2C-P3A (HAV)



---

Product Name	P2C-P3A (HAV)
Catalog No	IT-004-052Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in 10mM CBB, pH9.6, 0.1% SDS and 50% glycerol. Reconstitute the protein with 100 µl of Millipore water.
Protein	The E.Coli derived recombinant protein contains the P2C-P3A immunodominant regions, amino acids 1392-1521.
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.