


Human TIM-1 (Ectodomain)



Product Name	Human TIM-1 (Ectodomain)
Catalog No	IT-004-TIM-1p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of human TIM-1 (Ectodomain) protein(a.a. 17-294) (Genebank #: BAJ61033). TIM-1 is the member of T-cell immunoglobulin and mucin domain-containing proteins, which are the receptors for hepatitis A virus. hTIM-1 has recently been implicated as the potential cellular receptor for the Ebola virus.
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human TIM-1 (Ectodomain)