



# Anti-gD (HSV-1)

---

Product Name	Anti-gD (HSV-1)
Catalog No	IT-005-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Glycoprotein D of human herpes virus 1 (HSV-1) (Genebank No. AF372624)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HSV-1 glycoprotein gD. Cross-reactivity to gD from other herpes viruses not tested
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

