



# Anti-gB (HSV-2)

Product Name	Anti-gB (HSV-2)
Catalog No	IT-005-012M1
Source	Mouse monoclonal antibody, clone 8C9
Isotype	IgM
Formulation	Each vial contains 100 µg purified IgM in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing ectodomain of HSV-2 gB protein (Genebank No. ADG45140) (a.a. 23-767) (IT-005-012p).
Applications	ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with gB of HSV-2, but not HSV-1.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

