

Anti-gH/gL Complex(EBV/B95-8)



Product Name	Anti- gH/gL Complex(EBV/B95-8)
Catalog No	IT-005-030M8
Source	Mouse monoclonal antibody, clone 10B10
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Co-immunization with DNA vaccine expressing the ectodomains of gH protein (a.a.1-529) (Genebank #: P03231), and gL protein (1-137) (Genebank #: CAD53428).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with gH/gL complex(EBV/B95-8) and gH/gL/gp42 complex (EBV/B95-8)
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-005-030M8

