

gH (Ectodomain)(EBV)(Strain B95-8)



Product Name	gH (Ectodomain)(EBV)(Strain B95-8)
Catalog No	IT-005-030p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of gH (EBV) protein (a.a. 1-679) (Genebank #: P03231).
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gH (Ectodomain)(EBV)(Strain B95-8)