

Anti-NS1(Dengue Virus 1)



Product Name	Anti-NS1(DENV1)
Catalog No	IT-006-0050M2
Source	Mouse monoclonal antibody, clone31G9
Isotype	IgM
Formulation	Each vial contains 100 µg purified IgM in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing NS1(Dengue virus 1/GD05/99)(a.a. 1-353) (Genebank No. AAM12552)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Cross-reacts with NS1 from DENV2, 3, JEV, WNV, TEV.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-006-0050M2

