

NS1 (DV2)



Product Name	NS1 (DV2)
Catalog No	IT-006-0051p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged Dengue Virus-2 NS1 protein + V5 (Genebank No. FJ744720.1).
Applications	WB standard, antibody ELISA, vaccine, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NS1 (DV2) protein

Note: