

# Anti-NS1(TEV)



Product Name	Anti-NS1(TEV)
Catalog No	IT-006-0056M2
Source	Mouse monoclonal antibody, clone 29H3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing NS1(Tick-borne encephalitis virus) (a.a. 1-371) (Genebank No. NP_043135)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with NS1 from TEV. Not cross-reactive to NS1 from DENV1, 2, 3, 4, WNV, JEV and YFV as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-006-0056M2

