

NS1 (Tick-borne encephalitis virus)



Product Name	NS1 (Tick-borne encephalitis virus)
Catalog No	IT-006-0056p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged NS1 (Tick-borne encephalitis virus) protein (a.a.1-352) (GenBank Accession No. NP_775504).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

