



Anti-NS1 (Zika virus)

| | |
|----------------|--|
| Product Name | Anti-NS1 (Zika virus) |
| Catalog No | IT-006-0063m5 |
| Source | Mouse monoclonal antibody, clone 2C1 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol. |
| Immunogen | DNA vaccine expressing NS1 (Zika virus/Brazil)(IT-006-0063p). |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reactive to NS1 (Zika virus) only as tested. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

