



Anti-NS1 (Zika virus)

Product Name	Anti-NS1 (Zika virus)
Catalog No	IT-006-0063m8
Source	Mouse monoclonal antibody, clone 2F7A7
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing NS1 (Zika virus/Brazil)(IT-006-0063p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reactive to NS1 (Zika virus) and NS1 (St. Louis Encephalitis virus) only as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

