



# Anti-Human RBP-4

---

Product Name	Anti-Human RBP-4
Catalog No	IT-010-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Human RBP-4 protein (Genebank accession # <a href="#">NM006744</a> )
Applications	ELISA(1:10,000-1:500,000), etc.
Specificity	Reacts with Human ST2 protein.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.