




## P30 (SAG1) (Toxoplasma gondii)

Product Name	P30 (SAG1) (Toxoplasma gondii)
Catalog No	IT-011-001Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in 500mM NaCl, 50mM Phosphate buffer (pH 7.4), 200mM imidazole, 8M urea. Reconstitute the protein using 100 µl Ultrapure water.
Protein	C-terminal 6xHis tagged P30 (SAG1) immunodominant region (a.a. 45-198)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: P30 (SAG1) (a.a. 45-198)

Note: