

# Anti-VP1(Norovirus/Hu/GI.4/1643/2008/US)



Product Name	Anti-VP1(Norovirus/Hu/GI.4/1643/2008/US)
Catalog No	IT-013-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide.
Immunogen	DNA vaccine expressing VP1 of Norovirus(Hu/GI.4/1643/2008/US) protein(a.a.1-544) (Genebank accession # ACV41096).
Applications	ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with human Norovirus VP1 protein. Cross-reactivity to other family member not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

