

Anti-VP1 (Norwalk-like virus/ NLV/Saint Cloud/624/1998/US)



Product Name	Anti-VP1(Norwalk-like virus NLV/Saint Cloud/624/1998/US)
Catalog No	IT-013-003
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide.
Immunogen	DNA vaccine expressing VP1 of Norwalk-like virus (NLV/Saint Cloud/624/1998/US)) protein(a.a.1-540) (Genebank accession # AF414427).
Applications	ELISA(1:1000-1:100,000), WB, IP, Neutralizing assay, etc.
Specificity	Reacts with human Norovirus VP1 protein. Cross-reactivity to other family member not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-013-003

