

# Anti-GP(Marburg/Popp)



---

Product Name	Anti- GP(Marburg/Popp)
Catalog No	IT-014-002
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing ectodomain of GP(Marburg/Musoke) protein (a.a.24-649) (Genebank No. Z29337).
Applications	ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with most GP(Marburg) glycoproteins. Not cross-react with GP(Ebola).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

