

Anti-GP(Ebola/Cote d'Ivoire)



Product Name	Anti-GP(Ebola/Cote d'Ivoire)
Catalog No	IT-014-005
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	DNA vaccine expressing glycoprotein of (Ebola/ Cote d'Ivoire) (a.a.1-672) (Genebank No. NC_014372).
Applications	ELISA(1:1000-1:100,000), WB(1:200-1:2000), IP, etc.
Specificity	Reacts with GP(Ebola/Cote d'Ivoire) glycoprotein. Cross-reactivity with other subtypes not determined.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

