



# Anti-VP40 (Ebola/Zaire)

Product Name	Anti-VP40 (Ebola/Zaire)
Catalog No	IT-014-011M3
Source	Mouse monoclonal antibody, clone 5F3
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing VP40 (Ebola/Zaire) protein (GenBank No. Q77DJ6) (a.a. 1-326) (IT-014-011Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts only with VP40 (Ebola/Zaire) protein as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

