

# Anti-VP2(Canine parvovirus)



Product Name	Anti-VP2(Canine parvovirus)
Catalog No	IT-015-012EM1
Source	Mouse monoclonal antibody, clone 28F7
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant VP2 [Canine parvovirus] protein (a.a. 1-586) (Genebank Accession No. ACS50333)(IT-015-012Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Recognize VP2 [Canine parvovirus] protein. Cross-reactivity not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-015-012EM1**

