

Anti- Phlebovirus_G1 (Rift Valley Fever Virus)



Product Name	Anti- Phlebovirus_G1 (Rift Valley Fever Virus)
Catalog No	IT-017-001M1
Source	Mouse monoclonal antibody, clone 5G8
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Glycoprotein (Rift Valley Fever Virus) (IT-017-001Ep).
Applications	ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Phlebovirus_G1 (Rift Valley Fever Virus) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

