

# Anti-NP(SFTS Virus)(Huaiyangshan Virus)



Product Name	Anti-NP(SFTS Virus)(Huaiyangshan Virus)
Catalog No	IT-017-006M13
Source	Mouse monoclonal antibody, clone 14B2
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing NP (SFTS Virus)(Huaiyangshan Virus) protein (a.a. 1-245) (Genebank Accession No. AEA29891; or ADZ04472).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with NPs from SFTS Virus(Huaiyangshan Virus), RVFV, and Schmallenburg virus.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-017-006M13**

