




E2 (Rubella Virus)

Product Name	E2(Rubella Virus)
Catalog No	IT-018-003Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in 500mM NaCl, 50mM Phosphate buffer (pH 7.4), 200mM imidazole, 8M urea. Reconstitute the protein using 100 µl Ultrapure water.
Protein	C-terminal 6xHis tagged E2(Rubella Virus) immunodominant region (a.a. 31-105)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: E2(Rubella Virus) protein (a.a. 31-105)

Note: