

# Pr77 Gag Precursor (MMTV)



Product Name	Pr77 Gag Precursor (MMTV)
Catalog No	IT-019-001Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.
Protein	C-terminal 6xHis tagged mouse mammary tumor virus (MMTV) Pr77 Gag precursor protein (a.a. 1-591) (Genebank Accession No. AAA46543).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Pr77 Gag Precursor (MMTV)