


# Anti-Urease subunit gamma/beta [*Helicobacter pylori* J99]



---

Product Name	Anti-Urease subunit gamma/beta[ <i>Helicobacter pylori</i> J99]
Catalog No	IT-021-003
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Urease subunit alpha [ <i>Helicobacter pylori</i> J99] (Genebank accession # NP_222790)
Applications	Western Blotting (1:200-1:2000), ELISA, etc.
Specificity	Reacts with Urease subunit gamma/beta [ <i>Helicobacter pylori</i> J99]. Cross-reactivity to other subtypes not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Urease subunit gamma/beta [*Helicobacter pylori* J99]

Note: