


# Capsid (O'nyong-nyong virus)



Product Name	Capsid (O'nyong-nyong virus)
Catalog No	IT-022-005Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged O'nyong-nyong virus capsid protein (see below for its aa sequence, MW ~32kD).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



## Capsid Protein Sequence:

efipaqtyynrryqprpwtqrptiqvirpkpararaagqlaqlisavsrllartvpqkpararatkkqkqvkeqqstrnqkkk  
apkqkqtqkkkrpgrrrermcmkiendcifevkhegkitgyaclvgdkvmkpahvkgtidnadlaklafkrsskydleca  
qipvhmkdsaskfthekpegynwhhgavqysggfriptgagkpgdsgrpifdnkgrvvaivlgganegtrtalsvvtw  
nkdivtkitpegsve