


nsP3 (Ross River virus)



Product Name	nsP3 (Ross River virus)
Catalog No	IT-022-013Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	N-terminal GST and C-terminal 6xHis tagged Ross River virus nsP3 protein (a.a. 1333-1875)(Genbank Accession No. ACV66991).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>90% , based on SDS PAGE
Storage	Stable for 6 months at -20°C from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: nsP3 (Ross River virus), partial degradation is observed.