



Anti-N (HP-PRRSV)

Product Name	Anti-N (HP-PRRSV)
Catalog No	IT-023-001M3
Source	Mouse monoclonal antibody, clone 5A6
Isotype	IgG2a
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing nucleocapsid (HP-PRRSV) protein (IT-023-001Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with N (HP-PRRSV) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

