



## Anti-nsp7 (HP-PRRSV)

---

Product Name	Anti-nsp7 (HP-PRRSV)
Catalog No	IT-023-002M3
Source	Mouse monoclonal antibody, clone 3C9
Isotype	IgG2a
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing nsp7 (HP-PRRSV) protein (IT-023-002Ep).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with nsp7 (HP-PRRSV) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

