

Anti- envelope glycoprotein (Bourbon virus)



| | |
|----------------|--|
| Product Name | Anti- envelope glycoprotein (Bourbon virus) |
| Catalog No | IT-025-001M1 |
| Source | Mouse monoclonal antibody, clone 1C6 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol. |
| Immunogen | DNA vaccine expressing ectodomain of envelope glycoprotein (Bourbon virus) (IT-025-001p). |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with envelope glycoprotein (Bourbon virus) as tested. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

