

LecCapture™



Product Name	LecCapture™
Catalog No	IT-200-001
Source	Purified by affinity chromatography
Formulation	10 µg/ml Lectin in phosphate buffer, 100 ml in total.
Applications	Coating reagent for capturing glycoprotein for ELISA analysis. Use 100 µl/well for 96-well-dish.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.