


# Human HER3/ErbB3 Protein



Product Name	Human HER3/ErbB3 Protein
Catalog No	IT-300-005p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged Human c-erbB-3/c-erbB3/EEBB3/ErbB-3/erbB3-S/HER3/LCCS2/MDA-BF-1/p180-ErbB3/p45-sErbB3/p85-sErbB3 protein (a.a.20-643) (GenBank No. NP_001973).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Human HER3/ErbB3 Protein

Note: