

Anti- Human FGFR3 / CD333 (alpha(IIIb))



Product Name	Anti- Human FGFR3 / CD333 (alpha(IIIb))
Catalog No	IT-300-006M2
Source	Mouse monoclonal antibody, clone 4A11
Isotype	IgG2b
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing Human FGFR3 / CD333 Protein (alpha(IIIb)) Protein (IT-300-006p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Human FGFR3 / CD333 Protein (alpha(IIIb)) Protein as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

