


MOMP (C. trachomatis strain C/TW-3 Serovar C)



Product Name	MOMP (C. trachomatis strain C/TW-3 Serovar C)
Catalog No	IT-502-002Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea.
Protein	N-terminal 6xHis tagged MOMP (C. trachomatis strain C/TW-3 Serovar C) (a.a. 23-397) (GenBank Accession No. P08780).
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: MOMP (C. trachomatis strain C/TW-3 Serovar C)