


P1-IV adhesin (Mycoplasma pneumoniae M129)



| | |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | P1-IV adhesin (Mycoplasma pneumoniae M129) |
| Catalog No | IT-503-003Ep |
| Source | Recombinant protein purified from E.coli |
| Formulation | Each vial contains 100 µg protein in 50mM Tris-HCl (pH8.0), 150mM NaCl, 20% glycerol. |
| Protein | C-terminal 6xHis tagged P1-IV adhesin (Mycoplasma pneumoniae M129) (a.a. 1240-1522) (GenBank Accession No. NP_109829). |
| Applications | WB standard, antibody ELISA, immunogen, etc. |
| Purity | >95% , based on SDS PAGE |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: P1-IV adhesin (Mycoplasma pneumoniae M129)