



Anti-MB (Ureaplasma urealyticum)

Product Name	Anti-MB (Ureaplasma urealyticum)
Catalog No	IT-503-015M1
Source	Mouse monoclonal antibody, clone 2F12E5
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing MB (Ureaplasma urealyticum) (IT-503-015Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with MB (Ureaplasma urealyticum), cross-reacts weakly with MB (Ureaplasma parvum) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

