

# Human CCL5 / RANTES



---

Product Name	Human CCL5 / RANTES
Catalog No	IT-601-003Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in 50mM Tris-HCl and 150mM NaCl.
Protein	C-terminal 6xHis tagged Human CCL5 / RANTES (GenBank # NP_002976, aa 24-91)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

