



Anti-Human CD317

Product Name	Anti-Human CD317
Catalog No	IT-001-317
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Ectodomain of human CD317 (Tetherin) protein (a.a.24-162) (Genebank accession No.BC033873)
Applications	Western Blotting (1:200-1:2000), Flow cytometry
Specificity	Reacts with ectodomain of human CD317 (Tetherin) protein.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing human CD317

Note: