


# Human PPIB/Peptidyl-prolyl cis-trans isomerase B



Product Name	Human PPIB/Peptidyl-prolyl cis-trans isomerase B
Catalog No	IT-400-015p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human PPIB/Peptidyl-prolyl cis-trans isomerase B (aa 34-216; GenBank# P23284) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.
Reference	



E: Human PPIB

Note: