

Anti-Hemagglutinin Protein H (CDV)



Product Name	Anti-Hemagglutinin Protein H (CDV)
Catalog No	IT-007-003M6
Source	Mouse monoclonal antibody, clone 1A11C4
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Hemagglutinin Protein H (CDV) (IT-007-003p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with Hemagglutinin Protein H (CDV) as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

