

Anti-G [Nipah henipavirus]



Product Name	Anti- G [Nipah henipavirus]
Catalog No	IT-026-002M3
Source	Mouse monoclonal antibody, clone 6G2D12
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing glycoprotein [Nipah henipavirus] (IT-026-002p).
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with glycoprotein (Nipah henipavirus)(IT-026-002p), and weakly cross with glycoprotein (Hendra henipavirus)(IT-026-001p).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

