

# Human TrkA Protein/High affinity nerve growth factor receptor/NTRK1 isoform 4



Product Name	Human TrkA Protein/High affinity nerve growth factor receptor/NTRK1 isoform 4 / TrkA-III
Catalog No	IT-600-010p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	Human TrkA-III (aa 33 – 325; GenBank# P04629-4) with 6xHis tagged at C-terminus. This isoform does not bind NGF and does not interact with GRB2 and FRS2.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

