

Anti-Human ALK tyrosine kinase receptor



Product Name	Anti- Human ALK tyrosine kinase receptor
Catalog No	IT-400-030M4
Source	Mouse monoclonal antibody, clone 5B7F8
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Human ALK tyrosine kinase receptor.
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with Human ALK tyrosine kinase receptor as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

